



Respiratory ventilation products

Portable ventilators

LTV™ Series

LTV Series ventilators	4
LTV Series accessories.....	5
LTV Series circuits.....	6-10
LTV Series circuit accessories.....	11
LTV Series power accessories and batteries	12-15
LTV Series ventilator stands, packs and mounting accessories	16
LTV Series stands, packs and mounting accessories.....	17-19
LTV Series miscellaneous parts & accessories	20-21
LTV Series ventilator miscellaneous parts & accessories	22
LTV Series accessories, manuals & educational materials.....	23



**21102-001****LTV Series 1100 ventilator system**

LTV Series 1100 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter and operator manual

**18984-001****LTV Series 1150 system**

LTV Series 1150 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter, operator manual and quick reference guide

**18888-001****LTV Series 1200 system**

LTV Series 1200 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter, low flow oxygen adapter, operator manual and quick reference guide

**18888-201****LTV Series 1200 MR Conditional ventilator system**

LTV 1200 MR ventilator with protective boots, AC adapter with power cord, dust cover, circuit adapter, low flow oxygen adapter, LTV 1200 MR Conditional floor stand, case of 5 MR-safe 15 ft circuits and operator manual

**99503****LTV 1200 MR Conditional System with SprintPack™**

LTV 1200 MR Conditional ventilator, protective boots, AC adapter with power cord, SprintPack™ bracket, SprintPack battery power system, dust cover, circuit adapter, low flow oxygen adapter, LTV 1200 MR Conditional floor stand, case of 5 MR-safe 15 ft circuits and operator manual


LTM™ II Graphics Monitor and accessories



17650-001

LTM II Graphics Monitor package

LTM II Graphics Monitor, LTM mounting bracket, LTM II data cable, power splitter cable, compact flash card and LTM II operator manual



17916-001

Bracket, LTM Graphics Monitor mounting



11089

Cable, LTM data, 21 in



11089-002

Cable, LTM data, 6 ft



18093-001

Cable, power splitter, LTM II



18097-001

Stylus pen, LTM




18098-001

Card, compact flash

LTV Series 1100, 1150 and 1200 ventilator circuits without PEEP


Single patient use (SPU) circuits, 6 ft



29657-001

Patient circuit without PEEP, 22 mm SPU (10/pkg)


1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb



29658-001

Patient circuit without PEEP, 15 mm SPU (10/pkg)


1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb



29663-001

Patient circuit without PEEP with 1 water trap, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb with water trap, temp probe port, wye and 61 cm x 22 mm exp limb




29662-001

Patient circuit without PEEP with 1 water trap, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb with water trap, temp probe port, wye and 61 cm x 15 mm exp limb


Reusable circuits, 6 ft



17997-001

Patient circuit without PEEP, 22 mm reusable (ea) (replaces 10815)

1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb




18941-001

Patient circuit without PEEP, 15 mm reusable (ea)

1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb

LTV Series 1100, 1150 and 1200 ventilator circuits without PEEP (continued)


SPU heated wire circuits for Fisher & Paykel humidifiers, 5.5 ft



15090-102

Patient circuit without PEEP, 22 mm SPU single heated wire, Fisher & Paykel (10/pkg)


1.7 m x 22 mm heated insp limb, temp probe port, wye and 61 cm x 22 mm exp limb



15090-104

Patient circuit without PEEP, 15 mm SPU single heated wire, Fisher & Paykel (10/pkg)


1.7 m x 15 mm heated insp limb, temp probe port, wye and 61 cm x 15 mm exp limb



15091-102

Patient circuit without PEEP with elbow, 22 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 22 mm heated insp limb, temp probe port, wye, elbow and 61 cm x 22 mm heated exp limb




15091-104

Patient circuit without PEEP with elbow, 15 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye, elbow and 61 cm x 15 mm heated exp limb

SPU MR safe circuit, 15 ft




19189-001

Patient circuit without PEEP, 22 mm SPU, 15 ft (5/pkg)

LTV Series 900, 950 and 1000 ventilator circuits with PEEP


SPU circuits, 6 ft



29695-001

Patient circuit with PEEP, 22 mm SPU (10/pkg)


1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb



29710-001

Patient circuit with PEEP, 15 mm SPU (10/pkg)


1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb



29697-001

Patient circuit with PEEP with 1 water trap, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb with water trap, temp probe port, wye and 61 cm x 22 mm exp limb




29696-001

Patient circuit with PEEP with 1 water trap, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb with water trap, temp probe port, wye and 61 cm x 15 mm exp limb

LTV Series 900, 950 and 1000 ventilator circuits with PEEP (continued)


SPU heated wire circuits for Fisher & Paykel humidifiers, 5.5 ft



15090-101

Patient circuit with PEEP, 22 mm SPU single heated wire, Fisher & Paykel (10/pkg)


1.7 m x 22 mm heated insp limb, temp probe port, wye and 61 cm x 22 mm exp limb



15090-103

Patient circuit with PEEP, 15 mm SPU single heated wire, Fisher & Paykel (10/pkg)


1.7 m x 15 mm heated insp limb, temp probe port, wye and 61 cm x 15 mm exp limb



15091-101

Patient circuit with PEEP, 22 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 22 mm heated insp limb, temp probe port, wye and 61 cm x 22 mm heated exp limb




15091-103

Patient circuit with PEEP, 15 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye and 61 cm x 15 mm heated exp limb

LTV Series 800 ventilator circuits with PEEP and without flow sensor (for pressure triggering)


SPU circuits, 6 ft



29715-001

Patient circuit with PEEP, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb



29716-001

Patient circuit with PEEP, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb



10696
Adapter, circuit



17921-001
Exhalation valve with PEEP, reusable
(replaces 10824)



18172-001
Exhalation valve without PEEP, reusable
(replaces 10825)



17486-101
Diaphragm, exhalation valve
(for 17921-001 and 18172-001)

SprintPack™ Lithium-Ion Power System



19222-201

Power system, SprintPack lithium-ion with cable protector

SprintPack power manager, cable protector, 2 SprintPack batteries and operator manual



19333-201

Power manager, SprintPack with cable protector (Does not include SprintPack batteries)



19444-001

Battery, SprintPack



19096-001

Bracket, SprintPack mounting



19102-001

Transport pack for LTV and SprintPack
5 D-rings, detachable shoulder strap, hand strap and bedrail mount



11544

Cord, automotive power, 12 ft



19871-001

Mounting bracket, MR Conditional SprintPack



24555-001

SprintPack Cable Support Bracket

External sealed lead-acid batteries



10800

Battery kit, 9-hour external

12 V/33 amp-hour battery, battery box with strap, charger, battery cable and operator manual



10787

Battery, 9-hour external (12 V/33 amp-hour)



10790

Battery box with strap, external



10900

Battery kit, 3-hour external

12 V/12 amp-hour battery, padded battery bag with cable, charger, battery cable and operator manual



11475

Battery, 3-hour external (12 V/12 amp-hour)



10901

Cable, 3-hour external battery



11330

Bag with cable, padded battery (for use with 3-hour battery)



10802

Cable, DC external battery (for use with 3-hour and 9-hour battery)



10801

Charger, external battery (for use with 3-hour and 9-hour battery)

Universal power supply (UPS)



17894-001
Battery, universal power supply

Miscellaneous power accessories



11448
AC power adapter with power cord, LTV



17527-001
AC power adapter without power cord, LTV



10854
Fuse, 32 V, 10 amp battery



11456
Electrical strip with 4 m power cord,
4-outlet U.S. hospital grade



11544
Cord, automotive power, 12 ft



10536
Power cord, 110 VAC hospital grade, 10 ft

LTV Series stands



10611
Stand with oxygen cylinder ring,
LTV Series floor



10695
Stand, LTV Series table

LTV Series mounts, packs and carriers




10700
Mount, bedrail



11581-G
Backpack for LTV Series ventilators with
boots *(gray)*



19102-001
Transport pack for LTV Series ventilators
and SprintPack
5 D-rings, detachable shoulder strap, hand
strap and bedrail mount




19184-001
Freedom carrier, LTV Series ventilators and
SprintPack

LTV Series mounting accessories




10610
Bracket, female dove-tail



10979
Bracket, humidifier pole mount C-clamp



11454
Bracket, large battery mounting



17916-001
Bracket, LTM graphics monitor mounting



19096-001
Bracket, SprintPack mounting



10658
Hand strap, LTV Series ventilators
(gray)



19031-001
Hand strap, LTV Series ventilators
(purple)



11501
Support arm with circuit hanger, patient circuit

LTV Series mounting accessories (continued)



11503
Hanger, clip, circuit support arm



11450
Basket, utility



11451
Crossbar, short



11452
Crossbar, long



11453
Pole, humidifer



11455
Oxygen cylinder ring set
(for use with 10611)



17341-001
Strap, transport pack shoulder

LTV Series mounting accessories sold but not shown

Cat. no.	Description
CAS-1152	Caster for floor stand (for use with 10611)
CAS-1168	Caster for floor stand, locking (for use with 10611)

LTV Series mounting accessories *(continued)*



19870-001
L-bracket with screws for LTV Series MR Conditional floor stand



19871-001
Mounting bracket, LTV Series MR Conditional SprintPack



19873-001
Caster for LTV Series MR Conditional floor stand



19874-001
Caster for LTV Series MR Conditional floor stand, locking



19875-001
Oxygen cylinder ring set, LTV Series MR Conditional floor stand



34910
Cable, remote alarm, 50 ft
(for use with alarm 11553)



19103-001
Cable, remote alarm, BNC - phone, 20 ft



14647-001
Cable, nurse call splitter



10779
Cable, nurse call, normally closed



10780
Cable, nurse call, normally open



17986-001
Adapter, electrical, Fisher & Paykel 900 MR
900 dual heated limb




17987-001
Adapter, electrical, Fisher & Paykel 900 MR
901 single heated limb



10701
Adapter, low pressure oxygen



10702
Adapter, NIST oxygen



11550
Boot kit
(gray protective)



19175-001
Boot kit
(purple protective)



15016-001
Cap, blender



14446
Cap, tethered port
(for use with LTV Series 800, 900, 950, 1100
and 1150 ventilators)



18098-001
Card, compact flash



10109
Cover, front panel dust




11456
Electrical strip with 4 m power cord,
4-outlet U.S. hospital grade



10678
Fan grill assembly
Fan grill, fan grill enclosure, fan filter



10258x5
Filter, LTV Series air inlet (5/pkg)



10789x5
Filter, LTV Series fan (5/pkg)



10699
Hose, 15 ft oxygen
(green)



11466
Test lung, 1L



10934x50
Tie, circuit twist-lok, dual (50/pkg)



10936x50
Tie, circuit twist-lok, large (50/pkg)

LTV Series ventilator shipping carton



Cat. no.	Description
19818-001	Carton, shipping (<i>replaces 17496-001</i>)
17490-001	Sheet, divider (<i>secures AC adapter</i>)
17489-001	Sheet, divider (<i>secures operator manual</i>)
17487-001	Cap set, LTV end

Order all four parts

Cat. no.	Description
19802-001	Operator's Manual, LTV 1200/1150 English
19822-001	Quick Reference Guide, LTV 1200/1150 English
18999-001	In-Service DVD, LTV 1200 English
19990-001	Caregiver Training DVD, LTV 1150 English
19960-001	Emergency Set-Up DVD, LTV 1200 English
21106-001	Operator's Manual, LTV 1100 English
19795-001	Operator's Manual, LTV 1000/950/900 English
19819-001	Quick Reference Guide, LTV 1000/950/900 English
18581-001	In-Service CD/DVD, LTV 1000/950/900 English
18580-001	Caregiver Training DVD, LTV 1000/950/900 English
11600	Simulator CD, LTV 1000 English
17622-001	Recommended Settings Card, LTV 1000/950/900 English (25/pkg)
19797-001	Operator's Manual, LTV 800 English
19821-001	Quick Reference Guide, LTV 800 English
11524	Quick Reference Guide, HME/Hospital English
17754-001	Operator's Manual, LTM II Graphics Monitor English
11009	LTM Graphics Monitor CD, English
21188-001	Operator's Manual, LTV 1200/1150 Spanish
21196-001	Quick Reference Guide, LTV 1200/1150 Spanish
19799-001	Operator's Manual, LTV 1000/950/900 Spanish
19820-001	Quick Reference Guide, LTV 1000/950/900 Spanish