

# NewLife® Intensity® 10

High Flow Stationary Concentrator

CAIRE®  
*O<sub>2</sub> is what we do.*



**TIME-PROVEN DESIGN**  
**VERSATILE OXYGEN DELIVERY**  
**EASE OF SERVICE**



# Trusted therapeutic oxygen solution designed for medical equipment fleets world wide

## Built-to-last Oxygen Delivery Innovation

---

Since being introduced in 2009, the New Life Intensity 10 has been a trusted high flow stationary concentrator addressing patients' oxygen delivery needs in over 100 countries including specialty applications in home and medical care facilities.

### Time-proven design

- **Industry-leading reliability**—greater than 99% reliability rate out of the box and less than 2%\* warranty repair rate
- **Comprehensive warranty**—3-year warranty on the device
- **Durability**—engineered with an oversized compressor and sieve beds to minimize service

### Versatile oxygen delivery

- **Oxygen delivery**—up to 95.5% concentration at up to 10 LPM, meeting a wide range of patient prescriptions
- **High outlet pressure**—the first oxygen concentrator offering 20 psi to support specialty applications and up to 200 feet of oxygen tubing
- **Energy efficient**—reduced power consumption, minimizing the burden on the patient's energy bill

### Ease of service

- **Minimal maintenance**—proven quality design and performance to minimize required service
- **Modular design**—fewer components reduce time and effort to service the device
- **Enhanced diagnostics**—onboard alarm alerts enable quick troubleshooting

\*data on file



## NewLife Intensity 10 Specifications

### FLOW RATE

2.0 – 10.0 LPM

### DIMENSIONS

27.5 x 16.5 x 14.5 inches  
(69.9 x 41.9 x 36.8 cm)

### WEIGHT

58 lbs (26.3 kg)

### SOUND LEVEL

58 dB(A)

### OUTLET PRESSURE

20 psig (138 kPA)

### POWER

600 W

### MAXIMUM TUBING

200 ft (61 m)

### WARRANTY

3 years

### MAINTENANCE SCHEDULE

Felt Filter – 1 Year Replacement  
Intake Filter – Clean Weekly



## Accessories


### NASAL CANNULA

#### Adult

7 ft. (2.1 m) CU002-1  
25 ft. (7.6 m) CU002-4

#### Pediatric

7 ft. (2.1 m) CU002-7



### HUMIDIFIER BOTTLE

**Item Number** HU003-1

**Flow Rate 6-15 LPM** HU014-1



### MASK

#### Adult, Concentration

(Over the Ear Style)

7 ft. (2.1 m)  
CU006-1



#### Pediatric, Concentration

7 ft. (2.1 m)  
CU007-1



### OXYGEN TUBING

7 ft. (2.1 m) CU004-2  
25 ft. (7.6 m) CU004-3  
50 ft. (15.2 m) CU004-1

### OXYGEN OUTLET ADAPTOR

#### Item Numbers

**Black** F0025-1

**Green** B-776945-00



# SureFlow™

## Cost effective oxygen delivery for multiple patients from a single source

### Multi-patient Flow Control

Combine the NewLife Intensity 10 with the multi-patient SureFlow oxygen flow station to cost-effectively serve multiple patients with a single oxygen source. SureFlow is a uniquely engineered flow station that effectively administers medical-grade oxygen from an oxygen concentrator or another single oxygen source. SureFlow is configured with five low-flow flow meters, and can be used within a wide range of settings to accommodate a broad variety of patient needs, enabling oxygen delivery for critical care up to 10 LPM.

#### Positioning

- Set on a table or desktop
- Wall-mounted



#### Diverse applications

- Pediatric-to-adult use in hospitals, clinics, emergency shelters
- Underdeveloped areas with reduced to no infrastructure for therapeutic oxygen delivery
- Low flow pediatric oxygen needs in a critical care nursery in low to middle income countries

#### Ideal usage

- Serve up to five individuals per SureFlow unit
- Connect two or more SureFlow units together to serve more individuals

### SureFlow Specifications

#### FLOWMETERS

Low-flow flowmeters can be used from 1/10 LPM to 1 LPM with .1/10 LPM increments with back lines for proper viewing angle

#### FLOWMETER ACCURACY

±10% or 200mL whichever is greater

#### MAXIMUM INLET PRESSURE

20 psig (138 kPa)

#### DIMENSIONS

5.8 x 9.7 x 7.2 inches  
(14.7 x 24.6 x 18.3 cm)

#### WEIGHT

3.28 lb (1.49 kg)

#### WARRANTY

3 years on parts

#### INTERNATIONAL APPROVALS

CE 0459 (FM069-2 only)  
ISO 13485:2016-certified



Find out more about NewLife Intensity 10 and accessory products at

[CAIREinc.com](http://CAIREinc.com)

 CAIRE®

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2022 CAIRE Inc. All Rights Reserved. ML-CONC0084 H