

# HI FLOW STROLLER®

Liquid Oxygen Portables

## A safer way to breathe

Developed for oxygen patient transport when ambulation is needed, CAIRE, Inc.'s Hi Flow Stroller® is truly state-of-the-art cryogenics. This intrahospital system was designed to complement the existing hospital ambulatory oxygen source.

The Hi Flow Stroller is considerably easier to use and, in turn, recognized as safer for both the patient and hospital staff. As a practical, money-saving solution, the effectiveness and simplicity of the Hi Flow is the answer to portable oxygen needs.

A perfect choice for hospitals and patients requiring higher flow rates, the Hi Flow Stroller is also a natural fit for those patients who need a portable oxygen unit that can keep up with them. The higher oxygen flow capacity makes the Hi Flow Stroller a natural fit for activities like hiking and fishing.



## Features

- Ultra-lightweight: 8 lbs. (3.6 kg)
- 1.25 L capacity
- 20 psi operating pressure
- Available in comfortable black or grey soft case
- LED meter easy to read
- Side Fill and Top Fill configurations
- Increased ease of mobility



**CAIRE Contact Information:**  
Providers Toll Free  
1.800.482.2473

[www.caireinc.com](http://www.caireinc.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2019 CAIRE Inc. All Rights Reserved.

ML-LOX0015 F

# HI FLOW STROLLER®

Liquid Oxygen Portables

Product Specifications	
LOX capacity	3.0 lb (1,36 kg)
Gaseous equivalent capacity	1,025 L
Weight, filled	8.5 lb (3,86 kg)
Weight, empty	5.5 lb (2,49 kg)
Height	13.5" (343 mm)
Typical use time at 2 LPM	8 hours
Operating pressure	20 psi (137 kPa)
Normal evaporation rate	1.3 lb/day (0,57 kg/day)
Standard flow control range	Continuous Flow: Off, .5, 1, 2, 2.5, 3, 4, 6, 8, 10, 12, 15



## Ordering Information

Description	Catalog Number
Hi Flow SF (side fill) grey case	14089341
Hi Flow TF (top fill) grey case	14081315
Hi Flow SF (side fill) black case	14778770
Hi Flow TF (top fill) black case	14768299
Accessories	
Hi Flow Stroller Cart	20765565



**CAIRE Contact Information:**  
Providers Toll Free  
1.800.482.2473

[www.caireinc.com](http://www.caireinc.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2019 CAIRE Inc. All Rights Reserved.

ML- LOX0015 F